

BifurQFR

Bifurcation evaluation with Quantitative Flow Reserve and 3 Dimension Quantitative Coronary Angiography

Julien Adjedj, Shengxian Tu, Hans Reiber, William Wijns

Background

Virtual FFR

FFR

$$\text{FFR} = \frac{\cancel{(P_d - P_v)} / \cancel{R_s}^{\text{max}}}{\cancel{(P_a - P_v)} / \cancel{R_v}^{\text{max}}}$$

- Maximal hyperemia $\rightarrow R_s = R_v$
- $P_v \ll P_a$ et P_d

$$\text{FFR} = \frac{P_d}{P_a}$$

$$\text{FFR} = \frac{Q_s^{\text{max}}}{Q_N^{\text{max}}}$$

Virtual FFR/QFR



Algorithme

$$\frac{Dv}{Dt} = -\text{grad } V - \frac{1}{\rho} \text{grad } p + \frac{1}{\rho} \text{div } \tau$$

$$\tau_{ik} = \mu \left(\frac{\partial v_i}{\partial x_k} + \frac{\partial v_k}{\partial x_i} - \frac{2}{3} \delta_{ik} \frac{\partial v_s}{\partial x_s} \right) + \mu_\theta \delta_{ik} \frac{\partial v_s}{\partial x_s}$$






3D vessel modelling



Background


3D QCA





Tu et al. Int J Cardiovasc Imaging 2012.

Background

Bifurcation lesion



RAO 9, CRAN 63
DS 53%, 10.1 mm
Pro Bif Angle: 165
Dis Bif Angle: 48
Pro: 3.0 mm, 2.9 mm
Dis: 2.0 mm, 2.0 mm






QAngio XA 3D

Courtesy of Shengxian Tu

Background

QFR




Tu et al. JACC Cardiovasc int 2014

Tu et al. JACC Cardiovasc int 2016

BifurQFR

Objectives


- Use 3D QCA and FFR to evaluate predictors of SB suffering
- Evaluate QFR accuracy compared to FFR pre and post provisional T stenting



BifurQFR

Design


Observational & multicentric study





BifurQFR

Example

68 years old male
previous anterior STEMI



LAO 45, CAUD 51
DS 66%, LL 14.6 mm
Prox Bif Angle: 163
Dist Bif Angle: 35
PD: 2.8 mm, 2.7 mm
DD: 2.3 mm, 2.2 mm




QAngio XA 3D





BifurQFR

Example

68 years old male
previous anterior STEMI



LAO 32, CAUD 63
DS 42%, LL 7.3 mm
Prox Bif Angle: 170
Dist Bif Angle: 32
PD: 2.3 mm, 2.2 mm
DD: 2.3 mm, 2.1 mm





QAngio XA 3D

BifurQFR

Example



68 years old male
previous anterior STEMI



BifurQFR

Example

68 years old male
previous anterior STEMI








BifurQFR

Example

68 years old male
previous anterior STEMI



BifurQFR

Sample size and support

- **Sample size calculation: based on FAVOR trial**
- **90-100 patients with planned bifurcation PCI**
- **Distal LM are excluded**
- **Biotronik support using ORSIRO DES**
- **No FFR wire support**

Perspectives

- Better evaluate SB suffering predictors after provisional T stenting with 3D QCA
- Validate QFR in bifurcation lesions pre and post PCI



Thank you

julienadjedj@hotmail.com